

# Routine 3D shimming through automation

Problem: Regularly shimming magnets running in automation is time consuming

Solution: get iconNMR to do it for you

Implementation:

1. Copy your favourite shim file: 'wsh shim.icon'
2. In the iconNMR Lock/Shim Options -> Solvent/Probe Dependencies set the shim file for all solvents to 'shim.icon'
3. Download the au\_make\_shim.icon file from [http://nmr.chemistry.manchester.ac.uk/downloads/AU\\_programs](http://nmr.chemistry.manchester.ac.uk/downloads/AU_programs) and save it in /opt/topspin/exp/stan/nmr/au/src/user/. or equivalent
4. Create a standard  $^1\text{H}$  experiment and set the automation acquisition AU program 'AUNM' to au\_make\_shim.icon
5. Save the parameter set
6. Add the parameter set to your admin user in iconNMR
7. Put in a 90 %  $\text{H}_2\text{O}$ /10 %  $\text{D}_2\text{O}$  sample
8. Submit the experiment using iconNMR

What the program does:

1. Insert sample and lock
2. Wait 5 minutes for equilibration
3. Tune and match
4. 1D (z) gradient shimming
5. 3D gradient shimming
6. Save shim to shim.'todays\_date'
7. Save shim to shim.icon
8. Acquire a  $^1\text{H}$  spectrum

Takes 15-20 mins

Should run a separate line shape experiment afterwards

Generates a log of shims over time

# Routine 3D shimming through automation

```
/* au_make_shim.icon
```

```
Usage:
```

```
Use a 90:10 H2O:D2O sample
```

```
Equilibrates the sample then runs 3d shim and save a file to shim.icon and shim.todays_date
```

```
type XAUA au_make_shim.icon, or copy the 'PROTON' parameter set and replace "au_zg" with "au_make_shim.icon"  
*/
```

```
#include <stdio.h>
```

```
#include <time.h>
```

```
time_t rawtime;
```

```
struct tm *info;
```

```
char buffer[80];
```

```
time( &rawtime );
```

```
info = localtime( &rawtime );
```

```
strftime(buffer,80,"wsh shim.%Y%m%d", info);
```

```
/* equilibration time (seconds) five minutes should be enough*/
```

```
#define EQUILIBTIME 300
```

```
#define COMMAND1 ATMA
```

```
#define COMMAND2 XCMD("topshim tuneb tunea")
```

```
#define COMMAND3 XCMD("topshim 3d")
```

```
#define COMMAND4 XCMD(buffer)
```

```
#define COMMAND5 XCMD("wsh shim.icon")
```

```
#define COMMAND6 RGA
```

```
#define COMMAND7 ZG
```

```
/* Start of AU program */
```

```
GETCURDATA
```

```
i1 = ssleep (EQUILIBTIME);
```

```
COMMAND1
```

```
COMMAND2
```

```
COMMAND3
```

```
COMMAND4
```

```
COMMAND5
```

```
COMMAND6
```

```
COMMAND7
```

```
QUIT
```